


Product Datasheet

Golgi Antibody / GLG1 (orb1716243)

Description	The Golgi apparatus is an organelle present in all eukaryotic cells that forms a part of the endomem...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant full-length human protein was used as the immunogen for the Golgi Glycoprotein 1 antibody.
Storage	Aliquot the Golgi Glycoprotein 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the Golgi Glycoprotein 1 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q92896
Hazard Information	This Golgi Glycoprotein 1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (Frozen & FFPE): 0.5-1ug/ml (No special pretreatment is required for the immunohistochemical staining of formalin-fixed tissues)
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human breast carcin...