



Product Datasheet

GATA3 Antibody (orb1716253)


Description	GATA-3 (GATA binding protein 3) is a member of the GATA family of transcription factors. This 50kDa ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant full-length human protein was used as the immunogen for the recombinant GATA Binding Protein 3 antibody.
Storage	Aliquot the recombinant GATA Binding Protein 3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the recombinant GATA Binding Protein 3 antibody should be determined by the researcher.
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	P23771
Hazard Information	This recombinant GATA Binding Protein 3 antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human urothelial ca...



IHC staining of FFPE human lobular breas...



Western blot testing of human MCF-7 cell...