



Product Datasheet

GATA3 Antibody (orb1716255)


| | |
|----------------------------|---|
| Description | GATA-3 (GATA binding protein 3) is a member of the GATA family of transcription factors. This 50kDa ... |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Immunogen | Recombinant full-length human protein was used as the immunogen for the recombinant GATA Binding Protein 3 antibody. |
| Storage | Aliquot the recombinant GATA Binding Protein 3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles. |
| Note | For research use only |
| Application notes | Optimal dilution of the recombinant GATA Binding Protein 3 antibody should be determined by the researcher. |
| Formula | 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide |
| Isotype | Mouse IgG1, kappa |
| Clonality | Recombinant |
| Uniprot ID | P23771 |
| Hazard Information | This recombinant GATA Binding Protein 3 antibody is available for research use only. |
| Dilution Range | Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml |
| Expiration Date | 12 months from date of receipt. |



IHC staining of FFPE human urothelial ca...



IHC staining of FFPE human lobular breas...



Western blot testing of human MCF-7 cell...