



Product Datasheet

Oct4 Antibody / POU5F1 (orb1716300)

Description	OCT4/POU5F1 is a member of the POU family of transcription factors, involved in the regulation of pl...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant fragment from the human protein was used as the immunogen for the POU5F1 antibody.
Storage	Aliquot the POU5F1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Optimal dilution of the POU5F1 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q01860
Hazard Information	This POU5F1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human seminoma tiss...



IHC staining of FFPE human seminoma tiss...