

## Product Datasheet

### [KO Validated] TTC11/FIS1 polyclonal antibody (orb1716541)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | [KO Validated] TTC11/FIS1 polyclonal antibody  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Human, Mouse, Rat  |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | WB   |
| <b>Immunogen</b>           | Recombinant fusion protein of human TTC11/FIS1(NP_057152.2).                             |
| <b>Preservatives</b>       | 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2                              |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>MW</b>                  | 17KDa  |
| <b>Uniprot ID</b>          | <a href="#">Q9Y3D6</a>   |
| <b>Dilution Range</b>      | WB,1:500 - 1:2000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |