

Product Datasheet

Collagen I/COL1A1 polyclonal antibody (orb1717026)

Description	Collagen I/COL1A1 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	A synthetic peptide of human Collagen I/COL1A1(NP_000079.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	140KDa/220KDa
Uniprot ID	P02452
Dilution Range	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200 IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.