

Product Datasheet

Phospho-RPA2-S4/S8 polyclonal antibody (orb1717380)

Description	Phospho-RPA2-S4/S8 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S4 & S8 of human RPA2(NP_002937.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	36KDa
Uniprot ID	P15927
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.