

Product Datasheet

DiMethyl-Histone H3-K79 polyclonal antibody (orb1717522)

Description	DiMethyl-Histone H3-K79 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Immunogen	A synthetic dimethylated peptide around K79 of human histone H3(NP_003520.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	17KDa
Uniprot ID	Q16695
Dilution Range	WB,1:500 - 1:2000 ChIP,1:50 - 1:200
Expiration Date	12 months from date of receipt.