

Product Datasheet

NLRP5 polyclonal antibody (orb1717567)

Description	NLRP5 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human NLRP5(NP_703148.4).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	134kDa
Uniprot ID	P59047
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.