

Product Datasheet

MRP2/ABCC2 polyclonal antibody (orb1717792)

Description	MRP2/ABCC2 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human MRP2/ABCC2(NP_000383.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	200KDa-270KDa
Uniprot ID	Q92887
Dilution Range	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.