

## Product Datasheet

### SLC5A7 polyclonal antibody (orb1717828)

<b>Description</b>	SLC5A7 polyclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, WB
<b>Immunogen</b>	Recombinant fusion protein of human SLC5A7(NP_068587.1).
<b>Preservatives</b>	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	66-75KDa
<b>Uniprot ID</b>	<a href="#">Q9GZV3</a>
<b>Dilution Range</b>	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
<b>Expiration Date</b>	12 months from date of receipt.