

Product Datasheet

Caspr2/CNTNAP2 polyclonal antibody (orb1717996)

Description	Caspr2/CNTNAP2 polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Caspr2/Caspr2/CNTNAP2(NP_054860.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	160kDa
Uniprot ID	Q9UHC6
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.