

Product Datasheet

BTG3 polyclonal antibody (orb1718489)

Description	BTG3 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human BTG3.
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	29kDa
Uniprot ID	Q14201
Dilution Range	WB,1:200 - 1:1000 IHC,1:25 - 1:100
Expiration Date	12 months from date of receipt.