

Product Datasheet

DREB1B polyclonal antibody (orb1718577)

Description	DREB1B polyclonal antibody
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of arabidopsis thaliana DREB1B(NP_567721.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	25KDa
Uniprot ID	P93835
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.