

Product Datasheet

SARS-CoV-2 NSP10 polyclonal antibody (orb1718875)

Description	SARS-CoV-2 NSP10 polyclonal antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of coronavirus NSP10(YP_009725306.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	20-25KDa
Uniprot ID	PODTD1
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.