

Product Datasheet

ANAPC7 polyclonal antibody (orb1719113)

Description	ANAPC7 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human ANAPC7(NP_057322.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	67KDa
Uniprot ID	Q9UJX3
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.