

Product Datasheet

C1QTNF7 polyclonal antibody (orb1720253)

Description	C1QTNF7 polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human C1QTNF7(NP_001128642.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Refer to figures
Uniprot ID	Q9BXJ2
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.