

## Product Datasheet

### STAMBPL1 polyclonal antibody (orb1720344)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | STAMBPL1 polyclonal antibody   |
| <b>Species/Host</b>        | Rabbit   |
| <b>Reactivity</b>          | Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated   |
| <b>Tested Applications</b> | WB   |
| <b>Immunogen</b>           | Recombinant fusion protein of human STAMBPL1(NP_065850.1).                               |
| <b>Preservatives</b>       | 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2                              |
| <b>Concentration</b>       | 1mg/ml   |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Clonality</b>           | Polyclonal   |
| <b>MW</b>                  | 50kDa  |
| <b>Uniprot ID</b>          | <a href="#">Q96FJ0</a>   |
| <b>Dilution Range</b>      | WB,1:500 - 1:2000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.  |