

Product Datasheet

CEBP Delta/CEBPD polyclonal antibody (orb1720652)

Description	CEBP Delta/CEBPD polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human CEBP Delta/CEBPD(NP_005186.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	28KDa
Uniprot ID	P49716
Dilution Range	WB,1:200 - 1:2000
Expiration Date	12 months from date of receipt.