

Product Datasheet

MOCOS polyclonal antibody (orb1721139)

Description	MOCOS polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MOCOS(NP_060417.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	98kDa
Uniprot ID	Q96EN8
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.