

## Product Datasheet

# Peroxiredoxin 3 (PRDX3) polyclonal antibody (orb1721700)

<b>Description</b>	Peroxiredoxin 3 (PRDX3) polyclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, WB
<b>Immunogen</b>	Recombinant fusion protein of human Peroxiredoxin 3 (Peroxiredoxin 3 (PRDX3))(NP_006784.1).
<b>Preservatives</b>	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	28kDa
<b>Uniprot ID</b>	<b>P30048</b>
<b>Dilution Range</b>	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100
<b>Expiration Date</b>	12 months from date of receipt.