

Product Datasheet

TMEM126A polyclonal antibody (orb1721829)

Description	TMEM126A polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human TMEM126A(NP_115649.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	21kDa
Uniprot ID	Q9H061
Dilution Range	WB,1:500 - 1:2000
Expiration Date	12 months from date of receipt.