

Product Datasheet

MEP1B polyclonal antibody (orb1722502)

Description	MEP1B polyclonal antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human MEP1B(NP_005916.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	105kDa
Uniprot ID	Q16820
Dilution Range	WB,1:500 - 1:1000
Expiration Date	12 months from date of receipt.