

## Product Datasheet

### RPA70/RPA1 polyclonal antibody (orb1722547)

<b>Description</b>	RPA70/RPA1 polyclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IP, WB
<b>Immunogen</b>	Recombinant fusion protein of human RPA70/RPA70/RPA1(NP_002936.1).
<b>Preservatives</b>	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	68kDa
<b>Uniprot ID</b>	<a href="#">P27694</a>
<b>Dilution Range</b>	WB,1:500 - 1:2000 IP,1:50 - 1:100
<b>Expiration Date</b>	12 months from date of receipt.