

Product Datasheet

Methyl-HP1 gamma/CBX3-H175 polyclonal antibody (orb1722743)

Description	Methyl-HP1 gamma/CBX3-H175 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human HP1 gamma/CBX3(NP_009207.2).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Refer to figures
Uniprot ID	Q13185
Dilution Range	IHC,1:50 - 1:200
Expiration Date	12 months from date of receipt.