

Product Datasheet

MonoMethyl-Histone H2B-K5 polyclonal antibody (orb1722749)

Description	MonoMethyl-Histone H2B-K5 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	A synthetic monomethylated peptide around K5 of human H2BK5me1(NP_003519.1).
Preservatives	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Form/Appearance	1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	14kDa
Uniprot ID	Q16778
Dilution Range	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Expiration Date	12 months from date of receipt.