

## Product Datasheet

# DiMethyl-DNMT3A-K44 polyclonal antibody (orb1722750)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | DiMethyl-DNMT3A-K44 polyclonal antibody   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | IF  |
| <b>Immunogen</b>           | A synthetic dimethylated peptide around K44 of mouse DNMT3A(NP_031898.1).   |
| <b>Preservatives</b>       | 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2   |
| <b>Form/Appearance</b>     | 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2   |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Uniprot ID</b>          | <b>O88508</b>   |
| <b>Dilution Range</b>      | IF/ICC,1:50 - 1:200   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |