

## Product Datasheet

# Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) monoclonal antibody (orb1722779)

<b>Description</b>	Phospho-PAK4/5/6 (Ser474/Ser560/Ser602) monoclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A synthesized peptide derived from human Phospho-PAK4 + PAK5 + PAK6 (S474 + S560 + S602)
<b>Preservatives</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	Calculated MW: 64,80,75 kDa; Observed MW: 64,80,75 kDa
<b>Uniprot ID</b>	<a href="#">Q9P286</a> , <a href="#">O96013</a> , <a href="#">Q9NQU5</a>
<b>Dilution Range</b>	WB: 1/500-1/1000
<b>Expiration Date</b>	12 months from date of receipt.