

Product Datasheet

Acetyl-Histone H3 (Lys27) monoclonal antibody (orb1723265)

Description	Acetyl-Histone H3 (Lys27) monoclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IHC-Fr, IHC-P, IP, WB
Immunogen	A synthetic acetyl-peptide corresponding to residues surrounding Lys27 of human Histone H3
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 15 kDa; Observed MW: 15 kDa
Uniprot ID	P68431
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 ChIP: 1/20
Expiration Date	12 months from date of receipt.