

Product Datasheet

Phospho-p27 Kip 1 (Ser10) monoclonal antibody (orb1723760)

Description	Phospho-p27 Kip 1 (Ser10) monoclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser10 of human p27 KIP 1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 22 kDa; Observed MW: 27 kDa
Uniprot ID	P46527
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20
Expiration Date	12 months from date of receipt.