

## Product Datasheet

# Phospho-MEK1 (Thr292) monoclonal antibody (orb1723906)

<b>Description</b>	Phospho-MEK1 (Thr292) monoclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A synthetic phosphopeptide corresponding to residues surrounding Thr292 of human MEK1
<b>Preservatives</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	Calculated MW: 43 kDa; Observed MW: 43 kDa
<b>Uniprot ID</b>	<a href="#">Q02750</a>
<b>Dilution Range</b>	WB: 1/500-1/1000
<b>Expiration Date</b>	12 months from date of receipt.