

Product Datasheet

Phospho-MEK1 (Ser298) monoclonal antibody (orb1723907)

Description	Phospho-MEK1 (Ser298) monoclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser298 of human MEK1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Uniprot ID	Q02750
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200
Expiration Date	12 months from date of receipt.