

Product Datasheet

Phospho-FOXO3A (Ser253) monoclonal antibody (orb1724177)

Description	Phospho-FOXO3A (Ser253) monoclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser253 of human FOXO3A
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 71 kDa; Observed MW: 82 kDa
Uniprot ID	O43524
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.