

Product Datasheet

Phospho-CDK1/2/3 (Thr14) monoclonal antibody (orb1724396)

Description	Phospho-CDK1/2/3 (Thr14) monoclonal antibody
Species/Host	Rabbit
Reactivity	Hamster, Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Thr14 of human Cdk1
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Uniprot ID	P06493, P24941, Q00526
Dilution Range	WB: 1/500-1/1000 IP: 1/20
Expiration Date	12 months from date of receipt.