

---

## Product Datasheet

### Phospho-CDC6 (Ser106) monoclonal antibody (orb1724399)

**Description**

Phospho-CDC6 (Ser106) monoclonal

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	A synthetic phosphopeptide corresponding to residues surrounding Ser106 of human Cdc6
<b>Preservatives</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
<b>Form/Appearance</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>MW</b>	Calculated MW: 63 kDa; Observed MW: 63 kDa
<b>Uniprot ID</b>	<a href="#">Q99741</a>
<b>Dilution Range</b>	WB: 1/500-1/1000
<b>Expiration Date</b>	12 months from date of receipt.