

Product Datasheet

Phospho-Cannabinoid Receptor I (Ser316) monoclonal antibody (orb1724457)

Description	Phospho-Cannabinoid Receptor I (Ser316) monoclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser316 of human Cannabinoid Receptor I
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Uniprot ID	P21554
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.