

Product Datasheet

Phospho-Histone H4 (Ser1) polyclonal antibody (orb1724876)

Description	Phospho-Histone H4 (Ser1) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S1 of human Histone H4. Email for sequence
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 11 kDa; Observed MW: 17 kDa
Uniprot ID	P62805
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100
Expiration Date	12 months from date of receipt.