

---

## Product Datasheet

### Phospho-MSK1 (Ser376) polyclonal antibody (orb1725275)

**Description**

Phospho-MSK1 (Ser376) polyclonal

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC-P, IP, WB
<b>Immunogen</b>	A synthesized peptide derived from human Phospho-MSK1 (S376)
<b>Preservatives</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Form/Appearance</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 90 kDa; Observed MW: 90 kDa
<b>Uniprot ID</b>	<b>075582</b>
<b>Dilution Range</b>	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20
<b>Expiration Date</b>	12 months from date of receipt.