

Product Datasheet

Hsp90 alpha/beta polyclonal antibody (orb1725307)

Description	Hsp90 alpha/beta polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, IP, WB
Immunogen	A synthesized peptide derived from human Hsp90 alpha + beta
Preservatives	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 85 kDa; Observed MW: 90 kDa
Uniprot ID	P07900, P08238
Dilution Range	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100
Expiration Date	12 months from date of receipt.