

Product Datasheet

PSGR polyclonal antibody (orb1725455)

Description	PSGR polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human OR51E2. AA range:221-270
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Uniprot ID	Q9H255
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 ELISA: 1/10000
Expiration Date	12 months from date of receipt.