

Product Datasheet

Phospho-p107 (Thr369) polyclonal antibody (orb1725467)

Description	Phospho-p107 (Thr369) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	The antiserum was produced against synthesized peptide derived from human RBL1 around the phosphorylation site of Thr369. AA range:335-384
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P28749
Dilution Range	IHC: 1/50-1/100 ELISA: 1/10000
Expiration Date	12 months from date of receipt.