

Product Datasheet

Phospho-FGFR1 (Tyr654) polyclonal antibody (orb1725532)

Description	Phospho-FGFR1 (Tyr654) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human FGFR1 around the phosphorylation site of Tyr654. AA range:626-675
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 92 kDa; Observed MW: 120 kDa
Uniprot ID	P11362
Dilution Range	WB: 1/500-1/1000 IF: 1/50-1/200 ELISA: 1/10000
Expiration Date	12 months from date of receipt.