

## Product Datasheet

# Phospho-EDG1 (Thr236) polyclonal antibody (orb1725613)

<b>Description</b>	Phospho-EDG1 (Thr236) polyclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, WB
<b>Immunogen</b>	Synthetic peptide of human S1PR1
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 43 kDa; Observed MW: 42 kDa
<b>Uniprot ID</b>	<b>P21453</b>
<b>Dilution Range</b>	WB: 1/500-1/1000 ICC: 1/100-1/200 ELISA: 1/10000
<b>Expiration Date</b>	12 months from date of receipt.