

Product Datasheet

Phospho-DRP1 (Ser637) polyclonal antibody (orb1725629)

Description	Phospho-DRP1 (Ser637) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide of human DNM1L
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 82 kDa; Observed MW: 82 kDa
Uniprot ID	O00429
Dilution Range	WB: 1/500-1/1000
Expiration Date	12 months from date of receipt.