

Product Datasheet

Phospho-VEGF Receptor 2 (Tyr1175) polyclonal antibody (orb1725654)

Description	Phospho-VEGF Receptor 2 (Tyr1175) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	Synthetic peptide of human KDR
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P35968
Dilution Range	IHC: 1/50-1/100 ICC: 1/100-1/200
Expiration Date	12 months from date of receipt.