

Product Datasheet

Phospho-LAT (Tyr191) polyclonal antibody (orb1725678)

| | |
|----------------------------|---|
| Description | Phospho-LAT (Tyr191) polyclonal antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC-P, WB |
| Immunogen | The antiserum was produced against synthesized peptide derived from human LAT around the phosphorylation site of Tyr191. AA range:191-240 |
| Preservatives | Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3. |
| Concentration | 1mg/ml |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| MW | Calculated MW: 28 kDa; Observed MW: 28 kDa |
| Uniprot ID | O43561 |
| Dilution Range | WB: 1/500-1/1000 IHC: 1/50-1/100 ELISA: 1/10000 |
| Expiration Date | 12 months from date of receipt. |