

Product Datasheet

Phospho-GRLF1 (Tyr1087) polyclonal antibody (orb1725776)

Description	Phospho-GRLF1 (Tyr1087) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized phospho-peptide around the phosphorylation site of human GRF-1 (phospho Tyr1087)
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	Calculated MW: 171 kDa; Observed MW: 190 kDa
Uniprot ID	Q9NRY4
Dilution Range	WB: 1/500-1/1000 ELISA: 1/10000
Expiration Date	12 months from date of receipt.