

## Product Datasheet

### TLE1/2/3/4 polyclonal antibody (orb1725906)

<b>Description</b>	TLE1/2/3/4 polyclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P, WB
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from the C-terminal region of human TLE1/TLE2/TLE3/TLE4. AA range:721-770
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>MW</b>	Calculated MW: 83 kDa; Observed MW: 90 kDa
<b>Uniprot ID</b>	<a href="#">Q04727</a> , <a href="#">Q04725</a> , <a href="#">Q04724</a> , <a href="#">Q04726</a>
<b>Dilution Range</b>	WB: 1/500-1/1000 ELISA: 1/10000
<b>Expiration Date</b>	12 months from date of receipt.