

## Product Datasheet

# Hepatitis B Surface gen (HBsAg) (4H2) monoclonal antibody (orb1726142)

<b>Description</b>	Hepatitis B Surface gen (HBsAg) (4H2) monoclonal antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	HBsAg antibody was raised in mouse using highly purified hepatitis B surface antigen as the immunogen.
<b>Preservatives</b>	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P03141</a>
<b>Dilution Range</b>	ELISA: 1/10000
<b>Expiration Date</b>	12 months from date of receipt.