

Product Datasheet

Influenza A Hemagglutinin H7 (2H11) monoclonal antibody (orb1726145)

Description	Influenza A Hemagglutinin H7 (2H11) monoclonal antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Influenza A/Netherlands/219/03 (H7N7)
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.02% sodium azide, pH 7.3.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P03459
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.